Rev. A-02 2023-03-01

MINDFAB

MINDFAB GmbH Max-Josef-Metzger-Str. 6 86157 Augsburg, Germany
T.: +49 821 455252-0 | E-mail: augsburg@mindfab.com | www.mindfab.com

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FAB Cobalt-Chrome T4 CE 0425**LOT10000**

Cobalt based dental alloy for ceramic, type 4*

*According to ISO 22674:2016

The product is intended for dental laboratories for the creation of products that must be rigid and resistant, single crowns, bridges and superstructures on implants. The alloy FAB Cobalt-Chrome T4 is produced in conformity with standards ISO 9693:2020 and ISO 22674:2016, it is highly corrosion resistant, it does not contain Nickel and it is completely free of toxic elements like beryllium, indium and gallium. FAB Cobalt-Chrome T4 is a product intended to be used by qualified technician.

Intended purpose

Fabrication of metal base structures for metal-ceramic crowns and FPDs.

Intended patient population

This alloy is intended to be used for oral rehabilitation in edentulous or partially edentulous patients who suffer from loss of hard dental tissue. There are no limitations related to age, sex, clinical condition or other conditions, other than those possibly determined case by case by the professional user.

Intended users

The product is intended for professionals of the dental sector, in particular: prosthodontic laboratories, dental technicians, dentists. These professional figures shall have a high degree of specialisation and are perfectly aware of the destination of the product and the correct method of use for the realization of a final product and elaborate different methods to reach the best shape of the final prosthesis.

Contact and/or manipulation by other people is not envisaged.

Instruction for use**Modeling**

Modeling with a minimum thickness of 0.3mm on simple structures and 0.5mm in complex cases or with patients with bruxism. Avoid elbows.

Manufacturing

Separate the milled structures from the disk with tungsten carbide burs or with separating discs, clean with ultrasonic and acetone or bioalcohol to eliminate any oil residues from the working process.

Continue manufacturing process using tungsten carbide burs.

The worked framework has to be sandblasted with disposable equipment in aluminium oxide of 110 to 150 µm at a max. pressure of 3-4 bar and then steam clean. After cleaning, the framework should not be touched anymore and should be held with hemostatic forceps only.

Oxidation

Fire for 5 min. under vacuum at 950-980°C. After firing, the oxide layer has to be carefully sandblasted using disposable oxide-aluminium equipment 110 to 150 µm at a pressure of 2.5-4 bar. Steam and check that the surface has an homogeneous grey surface.

Re-use

The device can be reused until the available material is exhausted.

Ceramization

It's possible to use Bonding. Apply the ceramic, following its manufacturer's instructions. Slow cooling is recommended.

Soldering

We suggest using our MAGNUM SALDATURA Co to weld.

Waste disposal

The processing scrap must be disposed of as special waste in accordance with the EC directives 2008/98/CEE on waste, and 94/62/CEE on packaging and packaging waste and in compliance with national legislation in force on the subject.

Safety instructions

- Cobalt-Chrome-based alloys can seldom cause dermatitis on sensitive subjects. A Patch-Test is thus advisable.
 - Metal dusts and fumes are dangerous for health. Use an adequate exhaust fans during milling and polishing.
 - The alloy contains Cobalt, element classified as CMR (CAS Number: 7440-48-4).
 - Before prosthesis application verify if other metal implants are in patient's oral cavity. Coexistence of different metals can cause a "pile" effect.
 - We identify every batch with a number. We recommend to write it down in patient's file to allow its complete traceability.
 - This product does not need any special preservation precautions. MESA keeps batch traceability for 15 years.
 - The methods of preserving the products made with FAB Cobalt-Chrome T4 are under the responsibility of professional operators (dental technicians). This product must be stored, transported and used at a temperature between -30°C and +50°C to ensure the integrity of the packaging and the instructions contained therein.
 - MESA recommends that the patient be made aware of the possibility for dental alloys to affect MRI results.
 - The product is sold non-sterile.
- Users shall report any major accident occurred in relation to the use of this medical device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

A summary of safety and clinical performance (SSCP) is available on Eudamed.

Symbol	Description	Symbol	Description
	Batch number		Date of manufacture
	Consult instructions for use		
	Manufacturer		CE Marking
			Non-sterile
	Medical device		Catalogue number
	Temperature limitation		Contains hazardous substances
Percentage chemical composition (m/m)			
Co	66	Mo	6
Cr	27		
Others: Si, Mn			

FAB Cobalt-Chrome T4

Physical and mechanical properties	
Density	8.4 g/cm ³
Solidus-liquidus temperature	1307-1417°C
Melting point	1470°C
Yield load strength (Rp 0.2)	395 MPa
Percentage elongation at fracture	11%
Modulus of elasticity	233 GPa
Vickers hardness	255 HV10
Colour	White
Thermal expansion coefficient 25-500°C	14.3 · 10 ⁻⁶ K ⁻¹
Thermal expansion coefficient 25-600°C	14.5 · 10 ⁻⁶ K ⁻¹
Maximum cooking temperature	980°C
Suggested ceramics	VITA VM13
Cytotoxicity test according to ISO 10993-5	Passed