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## FAB Polyamide 12

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
### 2. GENERAL

FAB Polyamide 12 is a milling disc from microcrystalline synthetic (PA) material. Treatment of material is possible with all milling systems.

### 3. INDICATION

FAB Polyamide 12 is for use in the mouth of patient. It can be used as temporary construction as well as permanent solution or for orthodontic splints. Restorations are designed by a dental technician and produced by CAD/CAM technology.

### 4. CONTRAINDICATION

The product has no known contraindication.

### 5. SIDE EFFECTS

There is no known side effect of the product. There is no type IV irritation / allergy published.

### 6. PROCESSING

To process FAB Polyamide 12, use a soft material carbide cutter with a single cutting edge (USE PMMA MILLING BURS). FAB Polyamide 12, which has a different fracture toughness than PMMA, requires different milling parameters. Your work from FAB Polyamide 12 can be polished with a polishing paste and a soft brush. Please work only at low speeds to prevent overheating of the material during polishing.

### 7. MILLING

Use fine-toothed carbide burs or cutting discs to separate the milled restoration from the material. Cut the retaining bars between product and disc carefully and without applying pressure. In order to prevent plaque formation, a smooth polishing is mandatory. Polish with appropriate polishing pastes and polishing brushes.

## 8. CLEANING

After the preparation, the finished product is cleaned with compressed air. Steam cleaners are not recommended. Product deformation may occur at high temperature and pressure. Clean the product with water in an ultrasonic bath instead.

## 9. FIXING

Can be fixed with all temporary cements.

**ATTENTION: SPECIAL MILLING STRATEGY / PARAMETERS REQUIRED!!**

In this regard, please contact CAD/CAM manufacturer or MINDFAB GmbH.

## 10. PRODUCT DESCRIPTION

FAB Polyamide 12 is certified for use in the mouth of patient. They can be used as temporary construction as well as permanent solution or for orthodontic splints.

## 11. ADVANTAGES

- Very high material homogeneity thanks to the industrial polymerization process
- Well-balanced combination of fracture toughness and elasticity for clinical use. This guarantees excellent long-term stability
- Excellent abrasion resistance
- Ensures a high process safety, no mixing errors, no polymerization shrinkage and no unpleasant smell as no hand mixing or cartridge
- Economic fabrication due to digital workflow
- Lasting shades stability and aesthetics
- Natural fluorescence

## 12. TECHNICAL CUSTOMER SERVICE

### Manufacturer

on dent  2195  
 Karacaoglan Mah. Erena Plaza 6166 Sk., No: 32/1J  
 35070 Bornova / İZMİR, Turkey  
 T: 0 232 253 25 77 / 0 546 459 10 09  
 E-Mail: info@on-dent.com

### Global business partner

MINDFAB GmbH  
 Max-Josef-Metzger-Str. 6  
 86157 Augsburg, Germany  
 T: +49 821 455252-0  
 E-Mail: augsburg@mindfab.com

## 13. EXPLANATION OF SYMBOLS

  
 Consult instructions for use

 **LOT**  
 Batch Code

 **REF**  
 Catalogue number

  
 Manufacturer

  
 Keep away from sunlight / Keep in a dry place

  
 Use-by date

 **CE**  
 CE-mark

  
 Box contents

REF 000-1-001 REV 000 Print date: 02.05.2023