remanium[®]
star 

CE 0483

remanium[®] star powder 10 - 30 µm

remanium® star powder 10 - 30 µm

Pulver

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remanium® star powder 10 - 30 µm

Powder

CoCrW dental alloy based on cobalt for firing and partial denture technique, Type 5 acc. to DIN EN ISO 22674 / DIN EN ISO 9693 for processing with laser melting machines

1. Composition (% by mass)

Co	Cr	W	Si
60.5	28.0	9.0	1.5

Additional elements < 1%: Mn, N, Nb

This product contains cobalt.

Free from nickel, beryllium, lead, cadmium and iron.

2. Technical data

Yield strength*	$R_{p0.2}$	800 MPa
Tensile strength*	R_m	1170 MPa
Hardness*	H	395 HV10
Elongation at rupture*	A_5	11%; 13%**
Modulus of elasticity*	E	230 GPa
Density	ρ	8.6 g/cm ³
Solidus temperature	T_L	1320 °C / 2921 °F
Liquidus temperature	T_S	1420 °C / 3002 °F
CTE (25 °C - 500 °C / 77 °F - 932 °F)	α	14.4 x 10 ⁻⁶ K ⁻¹

The technical data represent typical values that are achieved when using remanium® star powder. A laser melting machine of the type Concept Laser Mlab cusing R 100 was used for the production of test specimens. The layer thickness was 25 µm.

* Variations of ± 10 % are possible depending on the batch.

** Typical values when using remanium® star powder with the sprint heat treatment.

3. Intended purpose

Non-precious alloy for the manufacture of medical products by means of an additive powder bed fusion process, e.g. selective laser melting.

4. Indication

remanium® powders for additive manufacturing are dental alloys based on cobalt. remanium® star powder is particularly suitable for the manufacture of fixed and removable prosthetic restorations, appliances and metal-ceramic frameworks.

5. Digital modeling

The dental products are modeled using a CAD software which has been approved for this purpose. Construction must be in line with dental regulations in order to fulfil mechanical requirements. The finished product must, however, have a material thickness of at least 0.3 mm.

6. Processing in the powder bed system

Before processing in the powder bed system begins, it must be clarified that the system is suitable for processing the material. There are product-specific, validated process parameters for the following machine manufacturers: Concept Laser GmbH, Trumpf Laser- und Systemtechnik

GmbH (Ditzingen, Germany), HBD (Guangdong Hanbang 3D Tech Co., Ltd.).

Further information can be found in our remanium® brochure „Processing instructions. Non-precious alloy for laser melting or milling“ (REF 989-814-00).

If you would like to work with a different supplier, please contact your machine manufacturer in order to clarify this or contact us. We will be pleased to help you check the system requirements.

7. Post treatment

Powder residues on the components manufactured by the additive manufacturing process must be removed with due care, taking safety precautions into consideration.

8. Stress relief heat treatment

The process outlined below describes stress relief heat treatment for products manufactured by additive manufacturing technology. The process ensures that the material properties are as indicated and guarantee that the components are of high quality and lasting durability. If the heat treatment program is altered, results may differ. This, in turn, may void the guarantee. Stress relief heat treatment optimizes the material structure of the CoCrW alloy. Inner stress is relieved and the fit of the dental component is improved. We recommend that heat treatment of the components is carried out under a protective atmosphere on a build plate in a suitable furnace. When using a gas box, observe the gas flow rates recommended by the manufacturer. Argon or nitrogen are suitable shielding gases. Heat treatment under ambient air causes a thicker reaction layer on the surface of the components which may negatively impact the optics and possibly restrict the function of the products.

Heat treatment classic:

1. Preheat the furnace to 600 °C / 1112 °F and then place the build plate with the supported components into the furnace.
2. Heat the furnace at maximum power to 1150 °C / 2102 °F.
3. Hold the temperature at 1150 °C / 2102 °F for 60 minutes, then switch off the heating.
4. Remove the build plate when the temperature has cooled down to 600 °C / 1112 °F.
5. Processing of the components can continue once they have cooled down.

Heat treatment speed:

1. Preheat the furnace to 1050 °C / 1922 °F and then place the build plate with the supported components into the furnace.
2. Hold the temperature at 1050 °C / 1922 °F for 60 minutes.
3. Then, remove the build plate and let it cool down at room temperature.
4. Processing of the components can continue once they have cooled down.

Heat treatment sprint (TruPrint 1000 and 2000):

1. Preheat the furnace to 1100 °C / 2012 °F and then place the build plate with the supported components into the furnace.
2. Hold the temperature at 1100 °C / 2012 °F for 45 minutes.
3. Then, remove the build plate and let it cool down at room temperature.
4. Processing of the components can continue once they have cooled down.

9. Post treatment

The components can be removed from the build plate using a separating disc or band saw. Support residues can be removed and the surface smoothed using a finely-toothed hard metal bur.

10. Powder processing

Excess powder can be passed through a sieve with a mesh of approx. 40 µm. The separated powder should undergo visual inspection and can, if considered in good order, be further processed in the powder bed system. Alternatively it can be stored, dry and under seal, in its container.

11. Welding

Suitable welding processes such as laser welding and TIG impulse welding produce a solder-free, mechanically strong and corrosion-resistant joint. It is important to observe the framework design, the surface structures, the welding sequence and the welding parameters recommended for each individual laser machine. Suitable welding wires are available as filler materials for all remanium® alloys. Examples: CoCr welding wire 0.25 mm (REF 528-215-10), 0.35 mm (REF 528-210-10) and 0.5 mm (REF 528-200-10).

12. Soldering

Soldering should be avoided if possible to reduce the number of materials. If, however, it is necessary to solder, please use a suitable solder and flux which are appropriate for the composition and melting range of the alloy used. Example: CoCrMo-Sold 1 (REF 102-306-00) in connection with rema® Flux 1 (REF 102-304-00).

13. Ceramic veneering

remanium® bonding alloys can be veneered with all suitable ceramic masses, e.g. ceraMotion® Me (Dentaurum), as long as the CTE is suitable for the alloy. The ceramic should cool down over a long period in the furnace unless the ceramic manufacturer has other specifications.

14. Contraindications and adverse reactions

Signs of intolerance to non-precious remanium® alloys are extremely rare if the manufacturer's Instructions for use are adhered to. If the patient has a proven allergy against any component within the alloy, this alloy must not be used for safety reasons. There have been individual reportings of local irritations which were electrochemically induced. If various alloys have been used, it is possible that galvanic effects may occur.

Exposure to dust and vapors may cause irritation to the eyes or respiratory tract. Individual alloy components may have a carcinogenic effect. Please adhere to warnings and safety data sheets. In the State of California elements of this alloy are considered to be carcinogenic.

15. Warnings and precautions

Metal dust may be produced as the alloy is mechanically processed. The mechanical process will also create heat and possibly burs. Processed material may therefore be hot and/or have sharp edges. It is possible that dust may be generated when the package is opened or the remanium® powder is transferred. Metal dust is combustible. Any type of ignition source should be kept at a distance. Metal fire powder extinguishers should always be positioned in close proximity to extinguish a fire if necessary. They should be in working order. Processing at temperatures above the solidus temperature may generate vapors. The workplace should be fitted with a suitable suction or ventilation system. We generally recommend that gloves, protective clothing, safety glasses and respiratory protection equipment are worn. Should there be contact with the skin, please wash the affected area with soap and water for a few minutes. Should the product come into contact with the eye, please rinse the eye

with a gentle stream of water or saline solution for at least 15 minutes. Please consult a doctor if you are exposed to individual elements.

Please observe the relevant safety data sheet (see www.dentaurum.com).

Interproximal or occlusal contact with other metallic elements and the use of different types of alloys in the same oral cavity should be avoided to prevent galvanic reactions or local irritations in individual cases. The alloy should not be ground or polished in the patient's mouth. There is no scientific evidence available on the safety or efficacy of treatment of pregnant women or nursing mothers or children.

16. Reuse

The alloy may be changed in its composition and properties if reused. In principle, we recommend using new material to manufacture high-quality prosthetic restorations.

17. Conditions for storing

Normal changes to the surrounding conditions (e.g. temperature, pressure, light) have no negative impact on the product characteristics of the remanium® alloys. remanium® powder should be stored in a dry place in its original packaging with the lid closed tightly.

18. Disposal

Please adhere strictly to the rules that apply for your region, bearing in mind the details outlined in the safety data sheets.

19. General information

Should the dental professional and/or the patient become aware of serious problems arising from the use of the product, it is important that

the manufacturer and the competent authority in the country in which the dental professional and/or the patient is resident are informed accordingly.

The SSCP is available at <https://ec.europa.eu/tools/eudamed> and www.dentaurum.com.

These Instructions for use contain the main steps for processing the non-precious remanium® alloys. Further information can be found in our remanium® star brochure „Processing instructions. Non-precious alloy for laser melting or milling“, REF 989-814-00.

For more information on our products and services, please visit www.dentaurum.com.

Our dental technical team is available to answer any questions that may arise. Dental Technical Hotline Tel No: + 49 72 31 / 803-410.

20. Explanation of symbols used on the label

Please refer to the label. Additional information can be found at www.dentaurum.com (Explanation of symbols REF 989-313-00).

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