GB Instructions for use dima Print Model

1. Introduction
2. Intended Use
3. Contra-indication
4. Hazard and Precaution
5. Storage conditions, expiry date and transport
6. Processing dima Print Model
7. Finishing
8. Plastic and packaging waste
9. Cleaning
10. Delivery units

The following instructions for use are intended for the user of dima Print Model as a material for dental situation models. This instruction for use provides also information about safety and environmental aspects; a safety datasheet is available on www.kulzer.com/downloads. In case more information is needed about the processing of dima Print Model material contact the Kulzer Service Hotline. Also see information at the end of this document. dima Print Model is a monomer based on acrylic esters for manufacturing of 3D-printed dental models. Suitable for printing situation models. dima Print Model shall be used in combination with cara Print 4.0.

2. Intended Use
dima Print Model is a monomer mixture based on acrylic esters for manufacturing situation models to be used in combination with cara Print 4.0 and HiLite power 3D.

3. Contra-indication
dima Print Model should not be used for any other purpose than situation models only. Any deviation from this instruction for use may have negative effect on the chemical and physical quality of dima Print Model. In case of an allergic reaction, please contact a medical physician.

4. Hazard & Precaution
Wear protection when handling dima Print Model. Protective glasses and nitrile gloves are advised. Information about the handling of the product can be found in the safety datasheet, which is available on www.kulzer.com/downloads.

5. Storage conditions, expiry date and transport
Store the product in the original packaging at room temperature in a dry and dark area, not exceeding 25°C (77°F). Close the packaging after each use. The expiry date of the product is mentioned on the product label. In case of exceeding the expiry date, the product is no longer guaranteed in terms of treatment. Do not expose to UV-light and moisture.

6. Processing dima Print Model
We advise to use nitrile gloves when using dima Print Model until post-curing.

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Before using dima Print Model, make sure to shake the product in the original packaging for approximately 5 minutes. Color deviations and print failures may occur when shaken insufficiently. Make sure that you work as clean as possible, dirty reservoirs or machines can cause deformation and therefore failure of the printed objects! Pour the liquid material in the reservoir of the 3D-printing machine. Follow the instructions for use for the 3D-printer to start the printer and the printing process. When the machine has finished its program remove the building platform from the machine as described in the instructions for use for the 3D-printer. Place the platform on some paper or cloth with the built jobs facing upwards. The printed jobs can now be removed from the platform using a plastic utility knife. Rinse the printed jobs twice with isopropanol to remove any excess material. Use an ultrasonic bath. Rinse the first time for 3 minutes, second time for 2 minutes. The second rinse must be carried out with clean isopropanol. Rinsing with alcohol solution must not take longer than 5 minutes, as this will cause defects in the printed job. After cleaning make sure the printed part is dry and free of solvent residues. Then place the printed jobs into the Kulzer HiLite power 3D for final polymerization. The final properties and also final color depend on the post-curing process. Post-curing is a blue-violet-light treatment to ensure that dima Print materials obtain full polymer conversion, through this the residual monomer content is reduced to a minimum and the highest mechanical properties are obtained. This procedure is a necessary step to produce a safe end-product. We strongly advice to use the Kulzer HiLite power 3D device.

Material	Total time (min.)	Post-curing device	Procedure
dima Print Model	6	HiLite power 3D	Turn objects after 3 minutes

The times for the post-curing and the turning of the object must be strictly followed. Please notice that the light sources of the light polymerisation device and the printing machine need routine maintenance. Follow the instructions for the relevant device.

7. Finishing
Remove any support structures and finish jobs if necessary, using conventional dental methods and instruments. Differences in color nuance may occur due to production in batches of the raw material and product, inadequate shaking of the original packaging before use or insufficient post-curing.

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Explanation of symbols on labelling	
	Batch code
	Manufacturer
	Keep away from sunlight
	Consult instructions for use
	Use-by date
	Storage temperature

8. Plastic and packaging waste
The product dima Print Model in its polymerized form is not harmful for the environment. Residual waste material in its liquid state should be delivered to a collection point for chemical waste material.

9. Cleaning
Polymerized dima 3D-printing materials should be cleaned with non-chemical products.

10. Delivery units
The product dima Print Model is available in the following packaging size: 1000 gr.

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