# Colado® CAD Ti2

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**Indikationen**
Colado CAD Ti2 ist suitable for permanent multi-unit dental restorations, e.g. bridges, according to DIN EN ISO 22674, type 3. In the anterior region, the bridge frameworks can be made with up to three consecutive pontics. The connector cross-section must not be less than 6 mm². In the posterior region, the bridge frameworks can be made with up to two consecutive pontics. The connector cross-section must not be less than 9 mm².

**Contraindikationen**
If the patient is known to be allergic or sensitive to titanium (grade 2), the material must not be used to fabricate restorations. Different types of alloy in the same oral environment can lead to a galvanic reaction.

**Processing Colado® CAD Ti2 after milling**
**1. Separating the frameworks from the disc**
The milled restorations may be separated from the milling disc using cross-cut tungsten carbide burs or separating discs suitable for processing titanium.

**Note:** Do not use any residue material from the discs as a casting alloy!

**2. Finishing / Cleaning**
Thoroughly clean the separated restorations with hot steam and subsequently degrease with ethyl alcohol in an ultrasonic bath. The restorations can be finished using clean tungsten carbide burs suitable for processing titanium so that no sharp angles or edges remain. Using the bur, grind the metal in one direction only in order to avoid overlapping and therefore prevent the likelihood of bubbles from occurring in the ceramic build-up. Please observe the manufacturer’s recommended maximum rpm for these instruments. Next, sandblast the surface using pure aluminium oxide (approx. 180 µm) at 2 – 3 bar pressure. Rinse off the jet medium under running water and then thoroughly clean the restoration with hot steam. After that, contamination of the object must be prevented (e.g. do not touch it with your fingers). When grinding the metal, do not inhale the grinding dust! Never use hydrofluoric acid!

**3. Veneering**
The metal wall thickness (single crowns 0.3 mm/ abutment crowns 0.5 mm) should always be observed. Colado CAD Ti2 frameworks can be veneered using the laboratory composite SR Nexco®/SR® Link (Ivoclar Vivadent) or a veneering ceramic suitable for titanium. Please observe the manufacturer’s instructions. The temperature during the ceramic firing must not exceed 800 °C (1472 °F), as this could cause a lattice transformation to take place within the metal structure.

**Caution:** Metal vapours and metal dust are harmful if inhaled; therefore, always use extraction equipment and/or suitable protective masks! When processing titanium and titanium alloys, there is a fire hazard. Therefore always use a cooling liquid when processing Colado CAD Ti2 in the milling machine.

**Storage**
No specific measures required.

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions. These regulations also apply if the materials are used in conjunction with products of other manufacturers.

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Dichte	4,51 g/cm³
Vickers hardness	150 HV 5/30
CTE value (25–500 °C)	9,6 x 10⁻⁶/K
Mass content of all components in mass <span> </span> %	Ti: ≥ 99% <p>Fe, C, O₂: ≤ 1.0%</p>
Tensile strength	≥ 345 MPa
0.2% proof stress	≥ 275 MPa
Fracture strain	20%
Modulus of elasticity	100 GPa
Melting point	1670 °C / 3038 °F

**Product Description**
Colado® CAD Ti2 are milling discs made from pure titanium (titanium grade 2) for the production of crowns and bridges using CAD/CAM technology. Colado CAD Ti2 fulfils the requirements according to DIN EN ISO 22674, type 3. The milling discs can be processed using conventional milling machines (please observe the manufacturer’s instructions).

**Note:** When using Colado CAD Ti2 8 mm, check the insertion holder (distance ring).

**Indications**
Colado CAD Ti2 is suitable for permanent multi-unit dental restorations, e.g. bridges, according to DIN EN ISO 22674, type 3. In the anterior region, the bridge frameworks can be made with up to three consecutive pontics. The connector cross-section must not be less than 6 mm². In the posterior region, the bridge frameworks can be made with up to two consecutive pontics. The connector cross-section must not be less than 9 mm².

**Contraindications**
If the patient is known to be allergic or sensitive to titanium (grade 2), the material must not be used to fabricate restorations. Different types of alloy in the same oral environment can lead to a galvanic reaction.

**Processing Colado® CAD Ti2 after milling**
**1. Separating the frameworks from the disc**
The milled restorations may be separated from the milling disc using cross-cut tungsten carbide burs or separating discs suitable for processing titanium.

**Note:** Do not use any residue material from the discs as a casting alloy!

**2. Finishing / Cleaning**
Thoroughly clean the separated restorations with hot steam and subsequently degrease with ethyl alcohol in an ultrasonic bath. The restorations can be finished using clean tungsten carbide burs suitable for processing titanium so that no sharp angles or edges remain. Using the bur, grind the metal in one direction only in order to avoid overlapping and therefore prevent the likelihood of bubbles from occurring in the ceramic build-up. Please observe the manufacturer’s recommended maximum rpm for these instruments. Next, sandblast the surface using pure aluminium oxide (approx. 180 µm) at 2 – 3 bar pressure. Rinse off the jet medium under running water and then thoroughly clean the restoration with hot steam. After that, contamination of the object must be prevented (e.g. do not touch it with your fingers). When grinding the metal, do not inhale the grinding dust! Never use hydrofluoric acid!

**3. Veneering**
The metal wall thickness (single crowns 0.3 mm/ abutment crowns 0.5 mm) should always be observed. Colado CAD Ti2 frameworks can be veneered using the laboratory composite SR Nexco®/SR® Link (Ivoclar Vivadent) or a veneering ceramic suitable for titanium. Please observe the manufacturer’s instructions. The temperature during the ceramic firing must not exceed 800 °C (1472 °F), as this could cause a lattice transformation to take place within the metal structure.

**Caution:** Metal vapours and metal dust are harmful if inhaled; therefore, always use extraction equipment and/or suitable protective masks! When processing titanium and titanium alloys, there is a fire hazard. Therefore always use a cooling liquid when processing Colado CAD Ti2 in the milling machine.

**Storage**
No specific measures required.

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions. These regulations also apply if the materials are used in conjunction with products of other manufacturers.

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