dima[®] Zirconia

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GB Instructions for use

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KULZER
MITSUBI CHEMICALS GROUP

Manufacturer:

Kulzer GmbH

Leipziger Straße 2

63450 Hanau (Germany)

CE 0197

dima[®] Zirconia ML Keramik-Rohlinge

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GB Instructions for use

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Indications

dima Zirconia discs are indicated for use as a substructure, for porcelain fused ceramic fixed dental restorations up to a maximum of 16 units. Limitations are listed in table 1.

Table 1: Indications for use and maximum number of pontics.

Process chain	Single unit crowns, inlays, onlays and veneers		Bridges (maximum number of pontics in between)*	
	Anterior	Posterior	Anterior	Posterior
dima Zirconia ML	X	X	2	2

* one free-end pontic up to premolar size is permitted

Contraindications

Contraindications are all applications not indicated in table 1 above.

Notes on construction

Wall and edge thickness	individual copings	bridges
Wall thickness	≥ 0,4 mm	≥ 0,5 mm
Edge thickness	≥ 0,2 mm	≥ 0,2 mm

Dimension requirements	anterior region	posterior region
Number of pontics	max. 2	max. 2
Connector cross section	6 mm²	9 mm²
Free-end pontic	–	max. 1
Connector cross section for cantilever	–	12 mm²

Application

- Fix the disc with the printed side (occlusal) upwards according to the CAM-System instructions, set the enlargement-factor.
- Start the machining.
- Do not use cooling lubricant or compressed air.
- After machining, dismount the blank with the framework.
- Remove the framework from the blank with caution using a diamond milling cutter.
- Gently remove dust from the framework, with a soft brush or an air-gun.*
- Position the dry framework in a sinter-bed inside the crucible.
- Programme your furnace, according to its manual and dima sinter instructions.
- Sinter the framework.
- Wait until cool down is complete before removing the sintered framework.
- Inspect the framework in respect of flaws, wall thickness, seat and margin fit. If necessary, make small adjustments e.g. using a water-cooled diamond-tool or another Zirconium dioxide suitable tool.
- Clean the framework with water and dry. It is now ready for veneering.

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GB Instructions for use

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		Ramp Rate °C/h	T _{max} °C	Dwell Time h	Cooling
	Crowns	480	1450	2	480°C/h
	Bridges	480	1450	2	Natural Furnace

Sinter instructions

		Ramp Rate °C/h	T _{max} °C	Dwell Time h	Cooling
dima Zirconia ML	Crowns	480	1450	2	480°C/h
	Bridges	480	1450	2	Natural Furnace

Safety instructions

- * Processing dima disc generates dust that can irritate the eyes, skin and airways. Therefore always ensure that the suction is working correctly.**
- Carefully read the MSDS.

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