

PMMA Blocks for Dental Use

Product overview





Table of contents

We are your partner	3
CADtools PMMA Clear	4 - 5
CADtools PMMA Mono	6 - 7
CADtools PMMA Multi	8 - 9
Shipping costs and payment options	10
Service and guarantee	11

All prices are valid per blank, plus statutory sales tax and shipping costs. Errors, misprints, price changes reserved.
Status: 05 / 2022



For your successful future

We are your partner

According to the vision as "leading counterpart for more added value in dental technology", we have laid a cornerstone, which directs actions in our company.

Best quality and a **fair price-performance ratio** combined with **transparency** and **simplicity** are our corporate principles. These values shape our strengths and create your success.

In our day-to-day cooperation, we concentrate on the fundamental issues that make an essential contribution to long-term value creation - economically and ecologically.

All these trend-setting guidelines serve one goal:
to make your work easier as a **long-term and reliable partner**, true to the Bavarian motto „ois isi“ (everything easy).

Just as we live and love our tradition in Bavaria, plastic materials have a corresponding tradition in dental technology. On the following pages, local delicacies show the properties of our materials



Transparent resin

CADtools PMMA Clear

With its maximum transparency, the PMMA Clear is ideally suited for all types of splints. The high level of transparency leaves space for the beauty of the teeth.

The CADtools PMMA Clear is as transparent as the Bavarian Purity Law for beer of 1516.

Chemical composition		Material data	
PMMA	99 %	Density	1,17 g / cm ³
Pigment	1 %	Flexural strength	> 120 MPa
		Vickers hardness	25 HV
		Water absorption	< 25 µg / mm
		Solution value	< 5 µg / mm ³

Dimensions	Article no. Clear	Price
98 x 16 mm	CLEAR-D98-16	16,50 €
98 x 20 mm	CLEAR-D98-20	17,00 €

✓ Your benefit:

- easy handling and processing
- reliable stability
- maximum transparency

→ Indications:

- up to 16 units
- provisional up to 36 months
- dental splint



Pre-colored resin

CADtools PMMA Mono

The pre-colored PMMA Mono saves the dental technician time and money and belongs in every dental laboratory.

Just like the gingerbread heart (Lebkuchen) with its colorful sugar decorations should not be missing at any Bavarian folk festival.

Chemical composition		Material data	
PMMA	99 %	Density	1,17 g / cm ³
Pigment	1 %	Flexural strength	> 120 MPa
		Vickers hardness	25 HV
		Water absorption	< 25 µg / mm
		Solution value	< 5 µg / mm ³

✓ Your benefit:

- easy handling and processing
- reliable stability
- demanding results through natural basic colors

→ Indications:

- up to 16 units
- provisional up to 36 months
- 16 color shades

Dimensions	Article no. A1	Article no. A2	Article no. A3	Article no. A3,5	Article no. A4	Article no. B1	Article no. B2	Article no. B3
98 x 16 mm	MONO-A1-D98-16	MONO-A2-D98-16	MONO-A3-D98-16	MONO-A3,5-D98-16	MONO-A4-D98-16	MONO-B1-D98-16	MONO-B2-D98-16	MONO-B3-D98-16
98 x 20 mm	MONO-A1-D98-20	MONO-A2-D98-20	MONO-A3-D98-20	MONO-A3,5-D98-20	MONO-A4-D98-20	MONO-B1-D98-20	MONO-B2-D98-20	MONO-B3-D98-20



Article no. B4	Article no. C1	Article no. C2	Article no. C3	Article no. C4	Article no. D2	Article no. D3	Article no. D4	Price
MONO-B4-D98-16	MONO-C1-D98-16	MONO-C2-D98-16	MONO-C3-D98-16	MONO-C4-D98-16	MONO-D2-D98-16	MONO-D3-D98-16	MONO-D4-D98-16	19,00 €
MONO-B4-D98-20	MONO-C1-D98-20	MONO-C2-D98-20	MONO-C3-D98-20	MONO-C4-D98-20	MONO-D2-D98-20	MONO-D3-D98-20	MONO-D4-D98-20	21,00 €

Pre-shaded multilayer resin

CADtools PMMA Multi

The all-rounder PMMA Multi gives even temporary solutions a high degree of aesthetics that is second to none. A real talent.

The pretzel is available in many variations: whether without soda wash, sweet or salty. Our PMMA Multi is just as versatile.

✓ Your benefit:

- easy handling and processing
- reliable stability
- demanding results through natural basic colors
- aesthetic color gradient

→ Indications:

- up to 16 units
- provisional up to 36 months
- 16 color shades
- natural color gradient

Chemical composition		Material data	
PMMA	99 %	Density	1,17 g / cm ³
Pigment	1 %	Flexural strength	> 120 MPa
		Vickers hardness	25 HV
		Water absorption	< 25 µg / mm
		Solution value	< 5 µg / mm ³

Dimensions	Article no. A1	Article no. A2	Article no. A3	Article no. A3,5	Article no. A4	Article no. B1	Article no. B2	Article no. B3
98 x 16 mm	MULTI-A1-D98-16	MULTI-A2-D98-16	MULTI-A3-D98-16	MULTI-A3,5-D98-16	MULTI-A4-D98-16	MULTI-B1-D98-16	MULTI-B2-D98-16	MULTI-B3-D98-16
98 x 20 mm	MULTI-A1-D98-20	MULTI-A2-D98-20	MULTI-A3-D98-20	MULTI-A3,5-D98-20	MULTI-A4-D98-20	MULTI-B1-D98-20	MULTI-B2-D98-20	MULTI-B3-D98-20



Article no. B4	Article no. C1	Article no. C2	Article no. C3	Article no. C4	Article no. D2	Article no. D3	Article no. D4	Price
MULTI-B4-D98-16	MULTI-C1-D98-16	MULTI-C2-D98-16	MULTI-C3-D98-16	MULTI-C4-D98-16	MULTI-D2-D98-16	MULTI-D3-D98-16	MULTI-D4-D98-16	49,00 €
MULTI-B4-D98-20	MULTI-C1-D98-20	MULTI-C2-D98-20	MULTI-C3-D98-20	MULTI-C4-D98-20	MULTI-D2-D98-20	MULTI-D3-D98-20	MULTI-D4-D98-20	54,00 €



Shipping costs and payment options

Shipping costs outside Germany may vary according to destination country

Price

UPS Express - End of Day 6:00 p.m.

Mo. - Fr.

9,49 €

Free shipping

From an order value of 500,- € (net) the delivery with UPS Express - End of Day 6:00 p.m. will be free of charge.

Payment options

Payment term (30 days net)

Our promise for your security

Service and guarantee



Satisfaction guarantee

Your 100 % satisfaction in the quality of the products is our promise. Should this not be the case, we will exchange the goods or refund the purchase price. At any time.



Refund warranty

We guarantee you fair and attractive conditions. If you find an identical product cheaper elsewhere within 6 weeks of the date of purchase, we will refund the difference.



3 year guarantee

We offer you extra security with our three-year guarantee on material properties and processing.



Overnight delivery

In case of placement of your order before 2 p.m., you will receive the goods the next working day (excluding peripheral areas).



Exchange service

You can exchange unopened goods or in mint condition for another product free of charge within 24 months.

from technician to technician

CADtools GmbH
Max-Josef-Metzger-Str. 6
86157 Augsburg | Germany

Phone: +49 821 455252-0
Mo. - Th. 8:00 am - 5:00 pm
Fr. 8:00 am - 4:00 pm

E-mail: augsburg@cadtools.eu



www.cadtools.eu

Businesspartner of CADdent®